

μCT at the Institute of Anatomy

David Haberthür¹, Ruslan Hlushchuk¹ and Valentin Djonov¹

¹Institute of Anatomy, University of Bern, Switzerland {haberthuer, hlushchuk, djonov}@ana.unibe.ch

{nabertnuer,nlusnchuk,ajonov}@ana.unlbe

Introduction

A the Institute of Anatomy of the University of Bern we have a μCT facility including two μCT machines; the SkyScan 1272 and 1172 (Bruker microCT, Kontich, Belgium, Fig. 1(b)).




Figure 1: SkyScan machines [1]



Figure 2: Form 2 3D printer [2]

This poster aims to introduce the μ CT facility at our institute and gives an overview over the capabilities and the scientific potential of the methodology. Additionally, we also have a high-quality, stereolithographic 3D printer (Formlabs Form 2, shown in Fig. 2) in our lab, which helps us to produce custom-made sample holders suited for scanning samples in wet and dry environment.

Materials and Methods

WE investigated different organs or entire lab animals, including zebrafish, mice, rats and mini pig.

Results






Figure 3: Visualization of a tomographic scan of a zebrafish, fixed in 4 % PFA.






(c) Inside view

Figure 4: Visualization of a μAngiofil-infused mini pig eye.

Results (continued)



(a) Cerebellum


1 mm

Aica

(b) Detail. ☆: Plexus choroideus

(c) Rectangle in subfigure (b)

Figure 5: Three-dimensional visualizations of μAngiofil-filled mouse brain.



(a) 3D Visualization

Figure 6: Visualization and quantification of number, size and distribution of lung metastases. Tomographic scan of an Durcupan-embedded mouse lung.

Discussion

CT is a powerful tool to gain insight into a broad kind of biological (and non-biological) samples. The new agent (μ Angiofil) shows excellent perfusion efficacy and continuous filling of the vascular and capillary bed. The μ angioCT-based morphometry protocols are superior to the time-consuming classical, histology-based morphometry. The high reproducibility combined with the aforementioned features makes μ angioCT with μ Angiofil a reliable partner in a very broad spectrum of basic research studies and clinically-oriented research.

Acknowledgments

THIS work has been supported by the Commission for Technology and Innovation grant 14055.1 (Dynamische Hochauflösende Mikroangiographie). We thank Dr. Helena Handrkova for providing the zebrafish and Adolfo Odriozola for the preparation of the lung sample.

References

- [1] Bruker microCT Products. URL: http: / bruker - microct . com / products / all _ products . htm (visited on 08/22/2016).
- [2] Formlabs Press Kit. URL: http://
 formlabs.com/company/press/
 (visited on 08/22/2016).